

Full Title: Updated Software for Automated assessment of glucose variability and quality of glycemic control in diabetes

Short running title: Automated calculations of glycemic variability metrics.

Authors:

V. Moscardó¹, M. Giménez², N. Oliver³, and N. R. Hill⁴.

Affiliation:

¹ Instituto Universitario de Automática e Informática Industrial, Universitat Politècnica de València, València, Spain

² Diabetes Unit, Endocrinology and Nutrition Department, Hospital Clínic Universitari, IDIBAPS, Barcelona, Spain

³ Faculty of Medicine, Imperial College, London, UK

⁴ Harris Manchester College, Mansfield Road, University of Oxford, United Kingdom

Word count: 2912

Number of tables and figures: 2

Corresponding author: Professor Nick Oliver

Division of Diabetes, Endocrinology & Metabolic Medicine

Imperial College London. St Mary's campus. Norfolk Place

London W2 1PG

Tel: + 44 (0)20 7594 2460

Fax: + 44 (0)20 7594 2432

Email address: nick.oliver@imperial.ac.uk

Abstract:

Background: Glycemic variability is an important factor to consider in diabetes management. It can be assessed with multiple glycemic variability metrics and quality of control indices based on Continuous Glucose Monitoring (CGM) recordings. For this, robust, repeatable calculation is important. A widely used tool for automated assessment is the EasyGV software. The aim of this work is to implement new methods of glycemic variability assessment in EasyGV and to validate implementation of each glucose metric in EasyGV against a reference implementation of the calculations.

Methods: Validation data used came from the JDRF CGM study. Validation of the implementation of metrics that are available in EasyGV software v9 was carried out and the following new methods were added and validated: PGS (personal Glycemic State), IGC (Index of Glycemic Control), times in ranges, and GVP (glycemic variability percentage). Reference values considered as gold standard calculations were derived from Matlab implementation of each metric.

Results: The Pearson correlation coefficient was above 0.98 for all metrics, except for Mean Amplitude of Glycemic Excursion (MAGE, $r=0.87$) as EasyGV implements a fuzzy logic approach to assessment of variability. Bland-Altman plots demonstrated validation of the new software.

Conclusions: The new freely available EasyGV software v10 (<https://www.phc.ox.ac.uk/research/technology-outputs/easygv>) is a validated, robust tool for analysing different glycemic variability and control metrics.

Key Words: Type 1 Diabetes, Glucose Variability, EasyGV software

Introduction

People with type 1 diabetes (T1D) face a life-long challenge to effectively manage glucose and minimise pre- and postprandial hyperglycemia while avoiding hypoglycemia (1). Glucose profiles in T1D can greatly differ even among people with HbA1c values approaching target, suggesting that glucose variability is an important component of dysglycemia (2). Reduction of glucose variability by addressing the risk of hypo- and hyperglycemia is a target in diabetes treatment, while variability has been associated with endothelial insult leading to macro- and microvascular complications (3,4,5,6), neuropathy (7), retinopathy (8), atherosclerosis (9), kidney disease (10), cardiovascular disease (11), and impairment of cognitive function (12).

Continuous glucose monitoring (CGM) provides immediate feedback on the glucose concentration, and on the magnitude, direction and rate of change of glucose in real time (13). These measurements enable direct assessment of the dynamics of glycemic fluctuations and calculation of variability metrics can give an objective reflection of homeostasis.

Numerous indices for the evaluation of glycemic variability are currently available, including CV (coefficient of variation), MAGE (Mean Amplitude of Glycemic Excursions), CONGA (Continuous overlapping net glycemic action), MAG (Mean Absolute Glucose), and the MODD (Mean of Daily Differences). Likewise, quality of glycemic control indices such as the M-value, GRADE (Glycemic Risk Assessment in Diabetes Equation), and the J-index; and the glycemic risk metrics such as LBGi (Low Blood Glucose Index), HBGi (High Blood Glucose Index) and the ADRR (average daily risk range) have also been described, though the J-index is a more robust measure of hyperglycemia by virtue of its tight correlations with HBGi, GRADE%_{hyper}, and the Hyperglycemia Index (14). Assessment of these measures has been implemented and made freely available in the EasyGV tool (15,16).

Newer metrics for the assessment of glycemic control based on CGM data have been defined for use in clinical practice and in the research setting. The Personal Glycemic State (PGS) is a composite index that assess four domains of glycemic control: mean glucose, glycemic variability, time in range and frequency and severity of hypoglycemia (17); the Index Glucose Control (ICG) transforms glucose values into a “score” with an adjustable

weight assigned to hypoglycemic values relative to hyperglycemic values (18); and the Glucose Variability percentage (GVP) has been defined (19).

All of these indices have been used in clinical research, but are not easy to derive and a validated platform is needed to make their calculation simpler, and to ensure homogeneity of calculation and reporting.

EasyGV (16) is a user-friendly software application developed in visual basic and enabled in an Excel (Microsoft) workbook which calculates glycemic variability metrics and quality of glucose control indicators from CGM data. This software platform was chosen as clinicians and researchers are familiar with it and it is intuitive, simple, and available across multiple platforms. Easy GV has 195 citations and 1719 downloads to date (correct in 2019).

Previously, EasyGV derived values for the followings metrics: MEAN, SD (Standard Deviation), CONGA, Lability Index (LI) , J-INDEX, LBGi, HBGI, GRADE, %GRADE-Hypo, %GRADE-hyper, %GRADE-Eu, MODD, MAGE, ADRR, M-VALUE, and MAG but had not been cross validated against a reference method . Moreover, the previous versions of EasyGV do not include the implementation of the newest above-mentioned metrics.

Several other glucose variability calculators are available – Glyculator (20), The MAGE computer program (21), The Glycemic Variability Analyzer Program (GVAP) (22), and CGManalysis (23) - that facilitate calculation of some, but not all, of the existing glycemic variability metrics. Table 1 summarizes their main features. EasyGV provides the most comprehensive and accessible assessment; we have therefore implemented the newer methods of glycemic variability assessment in EasyGV and have sought to validate implementation of each EasyGV metric against a reference implementation of the calculations.

Methods

a. EasyGV feature expansion

Three novel variability and quality of control glucose metrics - IGC (18), PGS (17), and GVP (19) - were included in the new version of EasyGV software v10. The mathematical formulae were taken from their original publications for inclusion in the updated version of the software. The following times in ranges have also been added: percentage of time below 50mg/dL, 54mg/dL, 70mg/dL and *Th mg/dL*; percentage of time in range 70 - 140

mg/dL , 70 - 180 mg/dL, and $Th1$ - $Th2$ mg/dL; and, percentage of time above 180mg/dL and Th ; where Th , $Th1$, and $Th2$ are customisable thresholds. A summary of the improvements carried out in EasyGV software v10 is included in the Supplementary material.

- MAGE calculation

The original description of the derivation of MAGE (24) has led to various interpretations of how to translate this into a computed algorithm. Indeed, there are several algorithms with small differences in their implementations (25,26). The MAGE calculation in EasyGV software v10 was carried out according to the original description as follows: Glycemic excursions exceeding 1 standard deviation (SD) were included and the direction of calculation was established by the direction of the first excursion (peak-to-nadir or nadir to-peak). Then, MAGE was implemented as this formula describes:

$$MAGE (mg/dL) = \frac{\sum AGE}{n}$$

where AGE is the amplitude of glycemic excursions and n is the number of glycemic excursions greater than 1 SD.

The MAGE calculation procedure is summarized in these main steps:

1. Determination of all local maximum and minimum values across the glucose recordings.
2. Evaluation of the maximum-minimum pairs values with the 1 SD criterion.
3. If the difference between the maximum-minimum pair is greater than 1 SD, the value is included in the global summation. If it does not, the pair difference value is disregarded.

b. Validation and data analysis

CGM data from 30 participants (10 adults, 10 adolescents, and 10 children) in the JDRF CGM study (<http://diabetes.jaeb.org/Dataset.aspx>) were selected at random to validate each variability and quality glucose control index in EasyGV v10. This sample size allows differences of 10mg/dL in the MAGE calculations to be identified between the reference standard and EasyGV v10 outcomes with 95% confidence and with 80% power.

Validation was carried out by comparison between the outcomes from Matlab implementation and the outputs of EasyGV v10. Matlab implementation of each metric

was based on the original publications; hence, the results obtained with this algorithm were considered the reference values for the EasyGV v10 validation.

The reference standard for MAGE calculation remains the manual 'pencil and ruler' approach. Therefore, validation was performed comparing MAGE manually calculated by an expert clinician (MG) with the outcomes of the MAGE-Matlab code. In this case, CGM data from 10 participants only was used with the duration of data for each participant at least 60 days.

The considered metrics in the updated EasyGV software v10, are: Mean; Standard Deviation (SD); CV; CONGA (27); LI (28); J-Index (29); LBGI and HBGI (30); GRADE (31) and the risk attributable to hypoglycemia and hyperglycemia from GRADE values (%GRADE-hypo, %GRADE-hyper) ; MODD (32); MAGE (24); ADDR (33); M-value (34); MAG (35); IGC (18); PGS (17); GVP (19); percentage of time between 70 and 140mg/dL, and between 70 and 180mg/dL; percentage of time below 50mg/dL, 54mg/dL, and 70mg/dL; and, percentage of time above 180mg/dL.

In EasyGV software v10 the user can change the interval used for CONGA; the default is 60 min or CONGA1, and this was used for the analysis. The LI period can also be changed; the default is 60min, and this was used in this analysis. For M-value calculations, the ideal glucose value can be set. The default of 120mg/dL was used here. In the IGC, the Upper and Lower limits of target range (ULTR and LLTR respectively), the exponents (a, b) and the scaling factors (c, d) are customisable. The default values were used in this analysis (ULTR= 80mg/dL; LLTR= 140mg/dL; a= 1.10; b= 2; c= 30; d= 30). In PGS, the upper and lower thresholds are configurable. For this analysis, 180mg/dL and 70mg/dL were used for the upper and lower threshold, respectively. Lastly, the user can calculate additional times above, below and within ranges.

Where there is missing CGM data, periods without glucose values that are longer than the defined 'Max Gap' are considered as gaps. Linear glucose interpolation can be included when the duration of missing data is less than the pre-defined Max Gap. This point is important when the CGM recordings are large since calibration periods or sensor changes could add error to the calculations.

c. Statistical analysis

Results were analysed by using Pearson Correlation in order to assess the consistency of the measurements made by both implementations (Matlab software and EasyGV v10). Furthermore, Bland-Altman plots for each metric were obtained to show the error of EasyGV v10 against the Matlab reference values. A p-value of <0.05 as considered to indicate statistical significance.

Results

a. Addition of new metrics

The implementation of PGS, IGC and GVP in addition to times in ranges and the hypoglycemia event detection function, which is needed for PGS calculations, was successfully undertaken. Figure 1 shows the instructions page (a) and main interface (b). The default value of the tuning parameters for each metric are seen. It is important to note that parameters a and b in IGC can only have values between 1 and 2.

In order to use EasyGV v10, CGM data are copied in to the Raw Data sheet. If glucose values are in mg/dL, the mg/dL >mmol/L from the “main” sheet has to be pressed before commencing the analysis. The *Max Gap* has to be defined. Depending on the defined max gap, the program will interpolate the data for each missing sample time, or it will analyse the raw data when the period without data is longer than the defined *max gap* time. The results are displayed in the “Results” sheet.

b. Validation

i. Manual MAGE vs MAGE-Matlab.

The analysis shows that there is consistency in the measurements made by both methodologies (Matlab and manually) measuring the same metrics ($r=0.9856$ with a 95% confident interval [0.9423; 0.9967] (p-value < 0.0001)) The mean difference is -0.0616 mmol/L and the Bland-Altman Plot for the mean differences shows randomness of the mean difference distribution.

Based on these results, the MAGE-Matlab implementation and the MAGE manual are equivalent. Therefore, the subsequent validation of EasyGV v10 with the Matlab implementation is feasible.

ii. EASY GV validation

All the correlations calculated for each metrics were 0.98 or better, except for MAGE where the Pearson correlation coefficient was 0.87. The statistical analysis is shown in the supplementary material with correlation plots between the reference and EasyGV metrics for all assessed glucose outcomes. The value of r approaching 1.0 across all evaluated metrics demonstrates that the implementation of Easy GV corresponds with the original and referenced definition of each one implemented in Matlab.

Discussion

The new metrics (PGS, IGC, GVP, and times in ranges) were successfully incorporated into EasyGV v10 and robust validation against reference, and in the case of MAGE, manual analysis.

The small discrepancies in the correlation analyses were not significant and may reflect the differing methods that each software uses to round the number of decimals considered to complete the calculations. For example, the discrepancies in GRADE may result from the decimal differences in the value of the double logarithmic function by both implementations. This error is then carried forward in the calculation of the several risk GRADE percentages. The same occurs with the PGS since this metric is obtained by means of the sum of four functions, and these functions have several further functions within their definition.

Of particular interest is the MAGE validation since, to date, the gold standard technique for this metric's estimation has been the manual pencil and ruler approach. Our software has been validated and has been shown to be comparable to the manual method. The small variance seen is due to the pre-processing of the temporal signal. The EasyGV v10 and Matlab software implementations follow slightly different techniques to filter the noise, further explaining the differences seen, but the nadir and peak detection was equivalent.

In addition, validation results for our MAGE implementation are similar to the results showed by other software in the literature (Fritzsche (21), GlyCulator (20), and Glycemic Variability Analyzer Program (GVAP) (22), which have been implemented to calculate MAGE. A comparative study (26) showed varying agreement among the available computer programs developed and validated to calculate MAGE. MAGE derived with EasyGV may differ from the values obtained with other calculators and it would therefore

be interesting to compare the results that EasyGV v10 provides with the results obtained with other software, along with the manual method, quantifying and identifying the sources of such differences where feasible. Of special interest is the comparison with Baghurst's algorithm (25), which has previously been assessed in a small dataset with short CGM periods and implemented in R resulting in a complete R package (23,25), but has not been validated against the pencil and ruler method.

Other calculators for a variety of variability metrics have been described. In addition to the variety of metrics considered, the main differences between the GlyCulator and our approach is the consideration of the missing data. EasyGV v10 considers a gap when there is lack of data during at least the minutes defined as a *max gap* by the user (50 mins is considered the default). However, the other software applies a linear interpolation. This approach may add some small errors in large CGM datasets. Fritzsche's MAGE computer programme has been assessed with short CGM recordings (72 hours) from a predominantly type 2 diabetes (T2D) cohort only. Lastly, both GlyCulator and GVAP consider fewer parameters and fewer that are customisable by the user than in EasyGV (e.g the *n* value for CONGA).

For optimal use of EasyGV, it is important to carefully check the data and correctly configure glucose units and the sampling frequency, as incorrect settings may give erroneous results. This is a weakness of the spreadsheet format and would be improved by implementing EasyGV 10 in solutions with direct access to raw glucose data, such as the Tidepool platform.

In addition, sampling times without glucose values should be entered in the spreadsheet as an empty cell to ensure robust results. The maximum duration of missing data for interpolation can be defined in the main spreadsheet of EasyGV (by default, it is defined as 50min). If gaps are greater than this threshold, data are considered internally as a new fragment of data and interpolation is not carried out. Not including time stamps in the spreadsheet is a further limitation of the EasyGV design but would increase complexity.

Another important point to consider is the upper and lower limit of CGM data, which are different depending on the different manufacturers. If glucose values are outside of the limit of detection, the sensor gives the saturation value (eg <2.2mmol/L). In these cases, values are considered constant at the saturated value (e.g. 2.2mmol/L) until the glucose

values return to the working region of the sensor. Different limits of detection may bias some reported values in EasyGV. There are several approaches to this issue. The simplest approach is to leave the data at the saturated value and report only those metrics unaffected by the potential differences, such as times above and below ranges which sit wholly within the continuous reporting range of the sensor. A second approach would be to replace values out of range with empty cells and report as missing data which will avoid biasing metrics of variability but may sacrifice the accuracy of times in ranges. Finally, a third approach would be to undertake a sensitivity analysis when analysing such data to assess the effect of either truncating the data or using the saturation point value.

Glucose variability assessment in T1D patients has become more relevant with the advent of continuous data and literature supporting its role in the pathogenesis of diabetes complications. The reduction of glucose variability is a potential therapeutic target and therefore robust, reproducible measurement is critical. The calculation of GV metrics is time consuming, has not been standardised, and is resource intensive. EasyGV v10 facilitates calculation to support glucose control evaluation and treatment. In the present work, a successful EasyGV v10 validation has been carried out as well as the incorporation of new indices and times in ranges. EasyGV v10 is, therefore, an appropriate tool to use for glycemic assessment in clinical practice and in clinical research. To our knowledge it is the only comprehensive, customisable and validated tool for glucose variability and continuous glucose data outcomes.

For MAGE, manual calculation is susceptible to human factor errors because the quality of calculations may be subject to the capacity and concentration of the person undertaking it. Thus, having an analysis tool that is able to undertake the operation with traceability to reference standard is critical.

EasyGV v10 computes all of the most important variability metrics and quality control indices in the current literature and can be customised, allowing the user to change the value of the tuning parameters.

One limitation of the EasyGV v10 is the fact that it needs a spreadsheet application or program to be used. However, this is unlikely to be an important limitation since most research groups and clinical teams have access to an office software package. Future work

might include implementation of the EasyGV software implementation in a web platform with the addition of graphical representations of reported metrics.

In summary, EasyGV v10 has been successfully validated to be a robust tool to use for glucose variability and control quality assessment. Moreover, the incorporation of novel metrics and the times in ranges in a simple interface make this tool useful, practical and comprehensive in glycemic evaluation.

Acknowledgements

V.M., N.R.H. and N. O directed the statistical analysis, interpreted the data, and wrote the manuscript. V. M. also carried out the validation process and the statistical analysis. M.G contributed to the process validation process, and interpreting the results and reviewing the manuscript. V. M and N. O. are the guarantors of this work and, as such, had full access to all of the data in the study and take the responsibility for the integrity of the data and the accuracy of the data analysis.

This paper presents independent research supported by the NIHR BRC at Imperial College Healthcare NHS Trust. The views expressed are those of the authors and not necessarily those of the NHS, the NIHR or the Department of Health. Vanessa Moscardó was recipient of a FPU grant, FPU13/04253, from the Spanish Ministry.

All the authors indicate there is no conflict of interest to disclose regarding this manuscript.

References:

- (1) Kovatchev, B. P. (2017). Metrics for glycaemic control—from HbA 1c to continuous glucose monitoring. *Nature Reviews Endocrinology*, 13(7), 425.
- (2) Monnier, L., Colette, C., & Owens, D. R. (2008). Glycemic variability: the third component of the dysglycemia in diabetes. Is it important? How to measure it?. *Journal of diabetes science and technology*, 2(6), 1094-1100.
- (3) Kilpatrick, E. S., Rigby, A. S., & Atkin, S. L. (2006). The effect of glucose variability on the risk of microvascular complications in type 1 diabetes. *Diabetes care*, 29(7), 1486-1490.
- (4) Hirsch, I. B., & Brownlee, M. (2007). The effect of glucose variability on the risk of microvascular complications in type 1 diabetes: response to Kilpatrick et al. and Bolli. *Diabetes care*, 30(1), 186-187.
- (5) Šoupal, J., Škrha Jr, J., Fajmon, M., Horová, E., Mráz, M., Škrha, J., & Prázný, M. (2014). Glycemic variability is higher in type 1 diabetes patients with microvascular complications irrespective of glycemic control. *Diabetes technology & therapeutics*, 16(4), 198-203.
- (6) Chen, X. M., Zhang, Y., Shen, X. P., Huang, Q., Ma, H., Huang, Y. L., ... & Wu, H. J. (2010). Correlation between glucose fluctuations and carotid intima-media thickness in type 2 diabetes. *Diabetes research and clinical practice*, 90(1), 95-99.
- (7) Bragd, J., Adamson, U., Bäcklund, L. B., Lins, P. E., Moberg, E., & Oskarsson, P. (2008). Can glycaemic variability, as calculated from blood glucose self-monitoring, predict the development of complications in type 1 diabetes over a decade?. *Diabetes & metabolism*, 34(6), 612-616.
- (8) Kilpatrick, E. S., Rigby, A. S., & Atkin, S. L. (2006). The effect of glucose variability on the risk of microvascular complications in type 1 diabetes. *Diabetes care*, 29(7), 1486-1490.
- (9) Snell-Bergeon, J. K., Roman, R., Rodbard, D., Garg, S., Maahs, D. M., Schauer, I. E., ... & Rewers, M. (2010). Glycaemic variability is associated with coronary artery calcium in men with Type 1 diabetes: the Coronary Artery Calcification in Type 1 Diabetes study. *Diabetic Medicine*, 27(12), 1436-1442.

(10) Văduva, C., Popa, S., Moța, M., & Moța, E. (2013). Analysis of Chronic Kidney Disease–Associated Glycemic Variability in Patients with Type 2 Diabetes Using Continuous Glucose Monitoring System. *Romanian Journal of Diabetes Nutrition and Metabolic Diseases*, 20(3), 315-322.

(11) Su, G., Mi, S., Tao, H., Li, Z., Yang, H., Zheng, H., ... & Ma, C. (2011). Association of glycemic variability and the presence and severity of coronary artery disease in patients with type 2 diabetes. *Cardiovascular Diabetology*, 10(1), 19.

(12) Rizzo, M. R., Marfella, R., Barbieri, M., Boccardi, V., Vestini, F., Lettieri, B., ... & Paolisso, G. (2010). Relationships between daily acute glucose fluctuations and cognitive performance among aged type 2 diabetic patients. *Diabetes care*, 33(10), 2169-2174.

(13) Kohnert, K. D., Heinke, P., Vogt, L., & Salzsieder, E. (2015). Utility of different glycemic control metrics for optimizing management of diabetes. *World journal of diabetes*, 6(1), 17.

(14) Rodbard, D. Metrics to Evaluate Quality of Glycemic Control: Comparison of Time in Target, Hypoglycemic, and Hyperglycemic Ranges with 'Risk Indices'. *Diabetes Technol. Ther.* **20**, 325–334 (2018)

(15) Hill, N. R., Oliver, N. S., Choudhary, P., Levy, J. C., Hindmarsh, P., & Matthews, D. R. (2011). Normal reference range for mean tissue glucose and glycemic variability derived from continuous glucose monitoring for subjects without diabetes in different ethnic groups. *Diabetes technology & therapeutics*, 13(9), 921-928.

(16) <https://www.phc.ox.ac.uk/research/technology-outputs/easygv>

(17) Hirsch, I. B., Balo, A. K., Sayer, K., Garcia, A., Buckingham, B. A., & Peyser, T. A. (2017). A simple composite metric for the assessment of glycemic status from continuous glucose monitoring data: implications for clinical practice and the artificial pancreas. *Diabetes technology & therapeutics*, 19(S3), S-38.

(18) Rodbard, D. (2009). New and improved methods to characterize glycemic variability using continuous glucose monitoring. *Diabetes technology & therapeutics*, 11(9), 551-565.

- (19) Garcia A, Balo AK, Buckingham BA, Hirsch IB, Peyser TA (2017). Application of Glycemic Variability Percentage: Implications for Continuous Glucose Monitor Utilization and Analysis of Artificial Pancreas Data. *Diabetes technology & therapeutics*, 19(2), 699-706
- (20) Czerwoniuk, D., Fendler, W., Walenciak, L., & Mlynarski, W. (2011). GlyCulator: a glycemic variability calculation tool for continuous glucose monitoring data. *Journal of diabetes science and technology*, 5(2), 447-451.
- (21) Fritzsche, G., Kohnert, K. D., Heinke, P., Vogt, L., & Salzsieder, E. (2011). The use of a computer program to calculate the mean amplitude of glycemic excursions. *Diabetes technology & therapeutics*, 13(3), 319-325.
- (22) Marics, G., Lendvai, Z., Lódi, C., Koncz, L., Zakariás, D., Schuster, G., ... & Tóth-Heyn, P. (2015). Evaluation of an open access software for calculating glucose variability parameters of a continuous glucose monitoring system applied at pediatric intensive care unit. *Biomedical engineering online*, 14(1), 37.
- (23) Vigers, T., Chan, C. L., Snell-Bergeon, J., Bjornstad, P., Zeitler, P. S., Forlenza, G., & Pyle, L. (2019). R Functions for Analysis of Continuous Glucose Monitor Data. *BioRxiv*, 625137.
- (24) Service, F. J. (2013). Glucose variability. *Diabetes*, 62(5), 1398.
- (25) Baghurst, P. A. (2011). Calculating the mean amplitude of glycemic excursion from continuous glucose monitoring data: an automated algorithm. *Diabetes technology & therapeutics*, 13(3), 296-302.
- (26) Sechterberger, M. K., Luijf, Y. M., & DeVries, J. H. (2014). Poor agreement of computerized calculators for mean amplitude of glycemic excursions. *Diabetes technology & therapeutics*, 16(2), 72-75.
- (27) McDonnell, C. M., Donath, S. M., Vidmar, S. I., Werther, G. A., & Cameron, F. J. (2005). A novel approach to continuous glucose analysis utilizing glycemic variation. *Diabetes technology & therapeutics*, 7(2), 253-263.

- (28) Ryan, E. A., Shandro, T., Green, K., Paty, B. W., Senior, P. A., Bigam, D., ... & Vantyghem, M. C. (2004). Assessment of the severity of hypoglycemia and glycemic lability in type 1 diabetic subjects undergoing islet transplantation. *Diabetes*, 53(4), 955-962.
- (29) Wójcicki, J. M. 'J'-index. A new proposition of the assessment of current glucose control in diabetic patients. *Horm. Metab. Res.* **27**, 41–2 (1995)
- (30) Kovatchev, B. P., Otto, E., Cox, D., Gonder-Frederick, L., & Clarke, W. (2006). Evaluation of a new measure of blood glucose variability in diabetes. *Diabetes care*, 29(11), 2433-2438.
- (31) Hill, N. R., Hindmarsh, P. C., Stevens, R. J., Stratton, I. M., Levy, J. C., & Matthews, D. R. (2007). A method for assessing quality of control from glucose profiles. *Diabetic medicine*, 24(7), 753-758.
- (32) Molnar, G. D., Taylor, W. F., & Ho, M. M. (1972). Day-to-day variation of continuously monitored glycaemia: A further measure of diabetic instability. Variation quotidienne de la glycémie enregistrée en continu: une mesure supplémentaire de l'instabilité diabétique. Tag-zu-Tag-Schwankungen des kontinuierlich gemessenen Blutzuckers: ein weiteres Maß für die Labilität des Diabetes. *Diabetologia*, 8(5), 342-348.
- (33) Kovatchev, B. P., Otto, E., Cox, D., Gonder-Frederick, L., & Clarke, W. (2006). Evaluation of a new measure of blood glucose variability in diabetes. *Diabetes care*, 29(11), 2433-2438.
- (34) Schlichtkrull, J., Munck, O. & Jersild, M. The M-Value, an Index of Blood-sugar Control in Diabetics (1965). *Acta Med. Scand.* **177**, 95–102.
- (35) Hermanides, J. *et al.* Glucose variability is associated with intensive care unit mortality. *Crit. Care Med.* **38**, 838–842 (2010)

TABLES

Table 1. Features of the main glycaemic variability calculation tools that are available in literature.

Glycaemic variability calculation tools	Features		
	Metrics available	Environment	Open source
EASYGV	<p>MEAN, Standard Deviation, CONGA, Lability Index, J-index, LBGI, HBGI, GRADE, %GRADE-Hypo, %GRADE-hyper, %GRADE-Eu, MODD, MAGE, ADRR, M-VALUE, and MAG.</p> <p>Times in ranges, PGS, ICG, GVP, have been included in the last version (EasyGV v10)</p>	User-friendly software application developed in visual basic and enabled in an Excel (Microsoft).	Yes
GLYCULATOR (20)	Standard deviation, %Coefficient of variation, MAGE, Weighted average of glucose values, J index, MODD, CONGA, Fractal dimension.	There are two available versions: Web based application and installable application.	Yes
THE MAGE COMPUTER PROGRAM (21)	MAGE	Software program	Yes

GVAP (22)	Average area above/below target range, Percentage spent above/below target range, CONGA, MODD, MAGE, Excursion Frequency, MAGE	Software application developed with Matlab	Yes
Cgmanalysis (23)	Sensor use, eA1c, GMI, median, quartiles, SD, CV, min/max, excursions above and below threshold, times in ranges, area under curve, MAGE, J-index, CONGAn, MODD, LBGI, HBGI	Implemented in R	Yes

Main Program Sheet

Options Sampling Interval 5 mins CONGA Length 60 mins LI Interval 60 mins M Value IGv 120 mg/dl MAGE CGM 1 (0=No; 1=Yes) Max Gap 50		PGS params Upper threshold 180 mg/dL Lower threshold 70 mg/dL IGC params LLTR 80 mg/dL ULTR 140 mg/dL a 1.10 range [1,2] b 2.00 range [1,2] c 30 d 30		hypoglycemia event definition duration 15 min separation 30 min %TIR - Range definition- Additional %TIR 0 (0=No; 1=Yes) Th1 0 mg/dL Th2 0 mg/dL %T below threshold - Threshold definition- Additional %T below Th 0 (0=No; 1=Yes) Th 0 mg/dL %T above threshold - Threshold definition- Additional %T below Th 0 (0=No; 1=Yes) Th 0 mg/dL	
--	--	--	--	--	--

Click to begin

(First step if data are in mg/dL)

1 of 1 Columns

EasyGV Version 10 Nathan R Hill - Copyright University of Oxford 2010-2016

Getting Started

Welcome to the EasyGV excel spreadsheet. This spreadsheet allows you to calculate 16 different measures of glycaemic variability (GV) using a simple interface. This workbook contains four sheets:

Instructions	This sheet detailing how to use the EasyGV calculator
Main	The sheet where you can set the various options of the GV methods and start the analysis
Raw Data	The sheet where you should input the data you want to analyse
Results	The sheet that will be produced when EasyGV has finished processing. A result of -999 means that method can not be calculated for the data

The excel workbook should work on Microsoft Excel versions 2007, 2010, 2011, 2013, 2016. It might work on other versions but they have not been tested.

The workbook uses macros to calculate the methods of GV so in order for it to work **you need to enable macros**.

Use the point as decimal separator (instead of comma).
EasyGV is free to use for academic and non-commercial purposes. If you wish to use it commercially please contact me on at nathan.hill@phc.ox.ac.uk
If you need technical support, you can also contact V. Moscardó at vmoscard@vamosgar@gmail.com

Options

The *Main sheet* contains a number of options their use is described here:

Sampling Interval	The time between each sample. The correct sampling interval will enable EasyGV to calculate profiles that are longer than 24 hours
CONGA Length	The CONGA length to use. See the CONGA paper to select an appropriate length the default is 60m
LI Interval	The Liability Index Interval to use default is 60m. See the Liability Index paper if you wish to change this
M Value IGV	The M-Value has an "index" of glycaemic variability the default is 120 but can be modified depending on which version of the M value you wish to use
MAGE CGM	The default is 0=off this means MAGE works as detailed in the service paper. BUT MAGE can work on CGM data Turn it on if you have CGM data
PGS parameters	The MAGE CGM contains a fuzzy logic algorithm that eliminates short term fluctuations due to sensor inaccuracy. Upper threshold (The default value is 180mg/dL). See the PGS paper if you wish to change this. Lower threshold (The default value is 70mg/dL). See the PGS paper if you wish to change this.
IGC parameters	ULTR (Upper Limit of Target Range, with default value of 140mg/dL). LLTR (Lower Limit of Target Range, with default value of 80mg/dL). a (exponent, in range from 1.0 to 2.0; with default value of 1.1). See the IGC paper if you wish to change this. b (exponent, in range from 1.0 to 2.0; with default value of 2.0). See the IGC paper if you wish to change this. c (scaling factor with default value of 30). See the IGC paper if you wish to change this. d (scaling factor with default value of 30). See the IGC paper if you wish to change this.
Hypoglycaemia event	duration (The default value is 15min)
definition	separation (The default value is 30min)
%TIR (customizable)	Th1 (lower threshold of the range)
%T below Th (customizable)	Th2 (upper threshold of the range)
%T above Th (customizable)	Th (customizable Threshold)
Interpolate	Th (customizable Threshold) You can interpolate any missing data you have using a straight line estimation. Maximum gap size is default at 50

Data Layout

The data should be laid out as follows. Missing data should just be left blank:

Header row	Patient 1
data	3.4
data	5
data	8
data	9
data	10
data	6
data	5.3
data	2.5
data	3.3
data	3.8
data	6
data	7.6

Glucose data should be in mmol/l. There is a convert to mmol/l button on the Main sheet to convert mg/dl to mmol/l

Figure 1. Appearance of the *Instructions* sheet (left) and the *Main* sheet (right).

Supplementary Material 1

Table 1. Correlation between Matlab and EasyGV software v10 outcomes for each assessed metric.

METRICS	r	CI (95%)	p-value
Mean	1.0000	[1.0000; 0.9999]	<0.0001
SD	1.0000	[1.0000; 1.0000]	<0.0001
CV	1.0000	[1.0000; 1.0000]	<0.0001
CONGA1	1.0000	[1.0000; 0.9999]	<0.0001
LI	1.0000	[1.0000; 0.9999]	<0.0001
JINDEX	1.0000	[1.0000; 1.0000]	<0.0001
LBGI	1.0000	[1.0000; 0.9999]	<0.0001
HBGI	1.0000	[1.0000; 1.0000]	<0.0001
GRADE	0.9874	[0.9740; 0.9939]	<0.0001
%GRADE-Eu	0.9998	[0.9996; 0.9999]	<0.0001
%GRADE-Hypo	0.9908	[0.9811; 0.9956]	<0.0001
%GRADE-Hyper	0.9998	[0.9997; 0.9999]	<0.0001
MODD	0.9982	[0.9962; 0.9991]	<0.0001
MAGE	0.8729	[0.7464; 0.9388]	<0.0001
ADDR	0.9980	[0.9959; 0.9990]	<0.0001
M ₁₂₀	1.0000	[1.0000; 1.0000]	<0.0001
MAG	0.9999	[0.9999; 0.9999]	<0.0001
IGC	0.9999	[0.9999; 0.9999]	<0.0001
PGS	0.9695	[0.9375; 0.9853]	<0.0001
GVP	0.9999	[0.9999; 0.9999]	<0.0001
%T above 180mg/dL	1.0000	[1.0000; 1.0000]	<0.0001
%T below 50mg/dL	1.0000	[1.0000; 1.0000]	<0.0001
%T below 54mg/dL	1.0000	[0.9999; 1.0000]	<0.0001
%T below 70mg/dL	1.0000	[0.9999; 1.0000]	<0.0001

%TIR (70-140mg/dL)	0.9991	[0.9982; 0.9996]	<0.0001
%TIR (70-180mg/dL)	0.9996	[0.9992; 0.9998]	<0.0001

Table 2. *Summary of the improvements carried out in EasyGV v.10 and its validation.*

METRICS	EasyGV software v. 9	EasyGV software v. 10
Mean	Correct metric implementation	Metric successfully validated
SD	Correct metric implementation	Metric successfully validated
CV	Not included	Metric included and successfully validated
CONGA1	Correct metric implementation	Metric successfully validated
LI	Correct metric implementation	Metric successfully validated
JINDEX	Correct metric implementation	Metric successfully validated
LBGI	Bug was identified	Corrected bug and metric successfully validated
HBGI	Bug was identified	Corrected bug and metric successfully validated
GRADE	Bug was identified	Corrected bug and metric successfully validated
%GRADE-Eu	Bug was identified	Corrected bug and metric successfully validated
%GRADE-Hypo	Bug was identified	Corrected bug and metric successfully validated
%GRADE-Hyper	Bug was identified	Corrected bug and metric successfully validated
MODD	Correct metric implementation	Metric successfully validated
MAGE	Correct metric implementation	Metric successfully validated

ADDR	Bug was identified	Corrected bug and metric successfully validated
M₁₂₀	Correct metric implementation	Metric successfully validated
MAG	Correct metric implementation	Metric successfully validated
IGC	Not included	Metric included and successfully validated
PGS	Not included	Metric included and successfully validated
GVP	Not included	Metric included and successfully validated
%T above 180mg/dL	Not included	Metric included and successfully validated
%T below 50mg/dL	Not included	Metric included and successfully validated
%T below 54mg/dL	Not included	Metric included and successfully validated
%T below 70mg/dL	Not included	Metric included and successfully validated
%TIR (70-140mg/dL)	Not included	Metric included and successfully validated
%TIR (70-180mg/dL)	Not included	Metric included and successfully validated

Figure 1 Representation of Matlab Outcomes versus EASYGV software v10 Outcomes for each assessed metrics.



